ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: Acute ankle sprains, chronic ankle instabilities, high ankle sprains, or prophylactic use. Ankle brace should be worn over an athletic sock. Extra Support (ES), Moderate Support (MS), and Light Support (LS) versions available.

CONTRAINDICATIONS: None.

WARNINGS AND PRECAUTIONS: If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately. Do not attempt to heat mold any component other than the foot plate. Do not use an oven, boiling water, or any heat source other than a heat gun. Refer to heat gun manufacturer's instructions for use for additional warnings and precautions. Do not heat foot plate to more than 300°F as melting may occur. If heated, the rivets will readily burn bare skin.

APPLICATION INFORMATION:

- A) Unfasten and loosen all straps. Place ankle in the brace and position at 90° (Fig. 1)
- B) Wrap and secure calf liner. (Fig. 2)
- C) Insert calf ratchet strap end through opening on the opposite calf ratchet strap. (Fig. 3) Pull, wrap, and lock both calf ratchet straps snugly into place without impairing circulation. (Fig. 4)
- D) Align tongue with the center of the forefoot. (Fig. 5, MS and ES versions only) Pull and wrap ankle strap(s) snugly into place without impairing circulation. (Fig. 6, MS and ES versions only)
- E) The lace guide on the tongue may be adjusted and the Achilles pad may be trimmed to fit individual variations in leg, ankle joint, and foot shape.

INSTRUCTIONS FOR HEAT MOLDABLE FOOT PLATE:

- A) Remove inner liner and outer shell (ES version) or ankle strap (MS version).
- B) Apply heat evenly to the area of interest using a heat gun. Do not attempt to heat mold any component other than the foot plate. Do not use an oven, boiling water, or any heat source other than a heat gun. Refer to heat gun manufacturer's instructions for use for additional warnings and precautions. Do not heat foot plate to more than 300°F as melting may occur. If heated, the rivets will readily burn bare skin.
- C) Manually form desired changes with gloved hands prior to cooling. Do not attempt to form components against bare skin.
- D) Reheat and continue forming as necessary.
- E) The foot plate will retain heat after heat source is removed. Allow sufficient cooling time prior to applying device.

CLEANING INSTRUCTIONS: To clean liners, hand wash in cold water using mild soap. Air dry. **Note:** If not rinsed thoroughly, residual soap may cause irritation and deteriorate material.

WARRANTY: Congratulations on your purchase of a DonJoy® brace. Please be assured that it is made of the finest material and technology available. Your Velocity™ Ankle brace comes with a 3 month warranty against defects on soft and rigid goods associated with the brace. Replacement soft goods can be purchased from our local representative, a medical provider, or by calling Customer Care.

Patent pending

For single patient use only

NOTICE: WHILE EVERY EFFORT HAS BEEN MADE IN STATE-OF-THE-ART TECHNIQUES TO OBTAIN THE MAXIMUM COMPATIBILITY OF FUNCTION, STRENGTH, DURABILITY AND COMFORT, THERE IS NO GUARANTEE THAT INJURY WILL BE PREVENTED THROUGH THE USE OF THIS PRODUCT.

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