



DONJOY®



THE CLEAR ADVANTAGE The IceMan® CLEAR³ cold therapy unit helps reduce pain and swelling, speeding up rehabilitation and recovery. The IceMan® helps provide extended cold therapy to aid a variety of indications and protocols directed by a medical professional. The IceMan® CLEAR³ utilizes DonJoy's patented recirculation technology, which helps maintain more consistent and accurate temperatures.



ICEMAN with dual pump recirculation technology. The DonJoy® IceMan® features a patented recirculation system that allows water warmed after flowing through the cold pad to be preserved and remixed with cooler ice water at a constant flow rate, providing consistent cool water distribution throughout the cold pad.



COMPETITORS Other versions restrict or pinch water flow to the cold pad, causing variable flow rates and inconsistent cold pad temperatures.

ICEMAN® CLEAR³

FEATURES

Proven Performance

- Time tested with over 20-years of experience delivering industry leading cold therapy solutions
- Developed through direct customer interaction and feedback

Technology Advantage

- Patented recirculation technology helps deliver consistent cold therapy throughout the cold pad, eliminating freezing inlet temperatures

Patient Convenience and Ease-Of-Use

- Cooler designed for easy set-up
- Plug-n-play operation
- Designed for easy portability and storage
- A complete line of cold pads to help deliver cold therapy to where it is needed the most

DONJOY® ICEMAN® CLEAR3

Description	Part Number
IceMan® CLEAR3 Cold Therapy Unit	11-0493
IceMan® CLEAR3 w/ Universal Cold Pad, S, RH	11-1629
IceMan® CLEAR3 w/ Universal Cold Pad, NS, RH	11-1630
IceMan® CLEAR3 w/ Universal Wrap-On Pad, NS, RH	11-1631
IceMan® CLEAR3 w/ Universal Wrap-On Pad, NS, EH	11-1633
IceMan® CLEAR3 w/ Shoulder, S Wrap-On Pad, NS, EH	11-1635
IceMan® CLEAR3 w/ McGuire Knee Wrap-On Pad, NS, RH	11-1636
IceMan® CLEAR3 w/ Ankle Wrap-On Pad, NS, RH	11-1638

WRAP-ON PADS

Description (Non-Sterile)	Part Number
Universal Wrap-On Pad, NS, RH (10.85" x 11.34")	11-0679-9-00000
Universal Wrap-On Pad, NS, EH (10.85" x 11.34")	11-0680-9-00000
Universal Wrap-On XL, NS, RH (12"x11.34")	11-0681-9-00000
Universal Wrap-On Pad XL, NS, EH (12"x11.34")	11-0682-9-00000
McGuire Knee Wrap-On Pad, NS, RH (10.85"x11.34")	11-1268-9-00000
Shoulder Wrap-On Pad S, NS, EH (12.75"x19"x12.25")	11-1397-9-00000
Shoulder Wrap-On Pad, NS, EH (14"x23"x12.5")	11-1334-9-00000
Ankle Wrap-On Pad, NS, RH	11-1522
Fore-Foot Wrap-On Pad, NS, RH	11-1521
Hip-Left Wrap-On Pad, NS, RH	11-9112
Hip-Right Wrap-On Pad, NS, RH	11-9111

COLD PADS

Description (Sterile)	Part Number
Universal Cold Pad, S, RH (10.85"x11.34")	11-0801-9-00000
Universal Cold Pad, S, EH (10.85"x11.34")	11-0802-9-00000
Universal Cold Pad XL, S, EH (12"x11.34")	11-0804-9-00000
Rectangle Cold Pad S, S, RH (4.81"x10.75")	11-0817-9-00000
Rectangle Cold Pad XL, L, S, EH (9.76"x13.75")	11-0813-9-00000
Wrist/Hand Cold Pad, S, RH (2.9"x4.2")	11-1269-9-00000
Ankle Cold Pad, S, RH (12.26"x8.28")	11-0815-9-00000
Description (Non-Sterile)	Part Number
Universal Cold Pad, NS, RH (10.85"x11.34")	11-0675-9-00000
Universal Cold Pad, NS, EH (10.85"x11.34")	11-0676-9-00000
Universal Cold Pad XL, NS, EH (12"x11.34")	11-0678-9-00000

STERILE DRESSINGS

Description (Non-Sterile)	Part Number
Universal Sterile Dressing	11-0822-9-00000
Universal XL Sterile Dressing	11-0823-9-00000
Rectangle XL Sterile Dressing	11-0820-9-00000
Shoulder XL Sterile Dressing	11-1360-9-00000

S = Sterile
NS = Non-Sterile
RH = Regular Hose (6.5")
EH = Extended Hose (15")

WARNING: THIS DEVICE CAN BE COLD ENOUGH TO CAUSE SERIOUS INJURY. DO NOT USE THIS DEVICE WITHOUT A PRESCRIPTION FROM A PHYSICIAN.



T 800.336.6569 D 760.727.1280 F 800.936.6569
 2900 Lake Vista Drive | Dallas, TX 75067 | U.S.A.
 DJOglobal.com/donjoy

Copyright © 2019 by DJO, LLC
 MKT00-8795 Rev A

Individual results may vary. Neither DJO Global, Inc. nor any of its subsidiaries dispense medical advice. The contents of this sheet do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.